VICAIR® ACTIVE **0**2

FOR THE ACTIVE WHEELCHAIR USER



Perfect for users who transfer frequently

Excellent lateral and frontal stability

Features a nifty storage pouch for personal items

The Active O2 offers the unique combination of ergonomically shaped foam and Vicair SmartCells. The four rear SmartCell compartments provide a high level of comfort and skin protection, and reduce pressure on the sitting bones and tailbone area. The foam front makes transfers easy and feels steady and improves wheelchair handling. The foam inlay is protected by an easy-toclean fluid-resistant liner. The rest of the cushion is 100% machine washable and breathable for optimal hygiene and perfect temperature & moisture regulation. The Active O2 has an easy accessible storage pouch for precious belongings (keys, phone, wallet). In short: everything you need to maintain your active lifestyle!

PRODUCT PERFORMANCE Skin protection / Anti decubitus A A A A A Pressure redistribution *** Shear force reduction *** Microclimate control Impact & vibration damping ____ Pocitioning

Positioning	
Lateral stability	4444
Frontal stability	4444
Sliding reduction	4444
Allows manual adjustments	4444
Allows easy transfers	4444

Maintenance & Safety	
Durability	AAAA
Reliability	AAAA
Functionality when misused	
Low maintenance level	AAAA
Level of hygiene	
Functionality & Application	***
Comfort	

4444

USED WITH

Spinal cord injuries, active users with a high transfer frequency

TECHNICAL SPECIFICATIONS	
Height	6 and 9 cm
Average weight	1250 grams
Max. user weight	200 kg
Machine washable	60°C (without foam)

CUSHION SIZES	
Standard	(Sizes >53 cm height 9 only)
On request	

